

Biopsy Needle Guides for Esaote



- Integral design without any accessories
- High grade medical stainless steel
- Artificial polishing
- Supporting accurate and angle adjustable biopsy
- Matching a wide range of needle gauges
- Simple pressing, quick release



About Us


As a China national high-tech enterprise, InnoFine are dedicated to R&D, manufacturing and sales ultrasound interventional products and consumables medical products since its establishment in 2013, providing safe and reliable, cost-effective products and overall solutions for thousands of hospitals in more than 100 countries around the world, products covering Ultrasound Needle Guides, Probe Cover, Gel, Biopsy Needle etc.

InnoFine will continue to adhere to the spirit of "Be innovative and forge ahead, Strive for excellence and perfection, Customer oriented ", and make unremitting efforts to promote the development of global interventional ultrasound field.

Cooperation Ultrasound Brand (Part)



Reference Sites (Part)



Classial Product (Part)



Biopsy needle guides for Esaote CA541 probe (Convex type probe)



Biopsy needle guides for Esaote P1-5 probe (Phased array type probe)



Biopsy needle guides for Esaote EC123 probe (Transvaginal type probe)



Biopsy needle guides for Esaote SC3123 probe (Microconvex type probe)



Biopsy needle guides for Esaote SL1543 probe (Linear type probe)



Biopsy needle guides for Esaote TRT33 probe (Biplane Transrectal type probe)

List of Biopsy Needle Guides for Esaote

REF#	Transducer	Tranducer Type	Gauge Size/ Guide Channel Length
JSM-002	CA421/CA430/CA431	Convex	14-22G,3-6Fr/3cm
JSM-106	CA621/CA631	Convex	14-22G,3-6Fr/3cm
JSM-053	CA541	Convex	11-23G,3-9Fr/3.5cm
JSM-198	AC2541/C1-8	Convex	11-23G,3-9Fr/3.5cm
JSM-232	CA711A	Convex	11-23G,3-9Fr/3.5cm
JSM-012	LA522/LA523/LA532	Linear	11-23G,3-9Fr/3.5cm
JSM-107	LA332/LA332E	Linear	11-23G,3-9Fr/3.5cm
JSM-197	SL1543/LA533/AL2443	Linear	11-23G,3-9Fr/3.5cm
JSM-253	AL2442/L3-11	Linear	11-23G,3-9Fr/3.5cm
JSM-113	LA435/LA424/SL2325	Linear	11-23G,3-9Fr/3.5cm
JSM-029	CA123	Microconvex	11-23G,3-9Fr/3.5cm
JSM-125	SC3123/CA1123/mC3-11	Microconvex	11-23G,3-9Fr/3.5cm
JSM-124	C5-2 R13/SC3121	Microconvex	18-20G/2cm
JSM-177	PA240	Phased Array	18-20G/2.5cm
JSM-192	PA230E	Phased Array	18-20G/2.5cm
JSM-196	SP2430/SP2730	Phased Array	18-20G/2.5cm
JSM-340	PA250/P1-5	Phased Array	14-22G,3-6Fr/3.5cm
JSM-193	TRT33/TLC3-13	Biplane Transrectal	14,16,18,20G/3cm
JEM-012	EC1123/SE3123(0°)	Transvaginal	16-18G/14cm
JEM-020	EC123/EC1123(3.8°)	Transvaginal	16-18G/13.5cm

(Please contact InnoFine for more models.)

Accurate, Safe and Efficient Operation



· Attach

· Angle Adjustment

· Needle Insertion

· Detach