



# SKE-134

ELECTRIC DIALYSIS CHAIR



# PRODUCT DATA SHEET



## SKE-134

### FEATURES AND HIGHLIGHTS



#### Technical parameters

▪ Size	1910×580mm
▪ Back	780×580mm
▪ Seat	480×580mm
▪ Leg	530×580mm
▪ Seat height	550-700mm
▪ Armrest	500×180mm
▪ Mattress thickness	70mm
▪ Back rest lifting	-8° - 70°
▪ Leg rest lifting	-50° - 0°
▪ Armrest swing out	0° - 270°
▪ Foot board adjustable distance	150mm
▪ Voltage and frequency:	100-240 V · 50/60 Hz
▪ Plug as required (hospital grade)	

#### Technical configuration

<input checked="" type="checkbox"/> With CPR	
<input checked="" type="checkbox"/> Hand controller	1pc
<input checked="" type="checkbox"/> Power cord	1pc
<input checked="" type="checkbox"/> 3" silent caster	1set
<input checked="" type="checkbox"/> Small dining table	1pc



# QUALITY CONTROL

SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating, efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality



**Steels**  
Baosteel from Fortune Global 500.



**Plastic**  
11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.




**Welding**  
Panasonic Robotic ensure 360° full smooth welding.




**Electric Motor**  
Linear motor, with UL/ROHS/EN standard. Noise under 50 dB, with lifetime 20,000 times movement.


## Accessories



**Chair Cover**  
Optional genuine leather or supebric cover, with artisanal making, comfortable hand feeling.



**Table Board**  
Adjustable ABS tray table, with antiskid convex edge, tank for cup and cellphone. USB/5V optional. Very convenient for donor.



**Handrail**  
Upgraded armrest, can be right-left, front-back adjustable. Donor can control universal position. Armrest is R120° concave curve, avoid donor's arm sliding, ensure the donor's safety and comfortable.



**Foot Board**  
150mm foot board adjustable distance.

## Ultrafiber leather color table



\* Note: the above is the physical scan, there may be a certain color difference in kind.



## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### South America Branch

**Saikang Medical Technology, Inc**  
W3400 Corrientes, Argentina  
TEL: +55 (11) 2758-0438  
Rua Universal, Hollywood,  
São Bernardo do Campo,  
Brazil  
TEL: +55 (11) 2758-0438

### East Africa Branch

**Saikang Medical Technology Ltd**  
Mezzanine Floor, Taj Tower,  
upper hill, nairobi, Kenya  
TEL: +254714084806

### South East Asia Branch

**Saikang Medical Technology Ltd**  
Sala Thammassop Subdistrict  
Thawi  
Watthana District  
Bangkok Thailand  
TEL: +66 642741046

### UK Branch

**Uk Healthward International Limited**  
DEPT 111,  
196 High Road Wood Green,  
London, N22 8HH  
United Kingdom  
TEL: +31 (0) 202111106

### Russia Branch

**LLC 'SAIKANG'**  
454004, Russian Federation,  
Chelyabinsk,  
Universitetskaya Nabe-  
reznaya  
TEL: +7 (919)300-00-06

### Mid East Branch

**Saikang Medical Supplies (Mid East) LLC**  
Musaffah  
P.O Box 53200, Abu Dhabi,  
UAE  
Tel: +971 585710652

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue