



# SKH046-11

OVERBED TABLE




**SAIKANG** / Medical

Healthier  
your life

# SKH046-11

## FEATURES AND HIGHLIGHTS



### Technical parameters

External size (LxWxH) 800×395×870-1170mm

### Technical configuration

- Air Spring for Height Adjustment 1pc
- Phenolic Resin Table Top 1pc

#### Material

Powder coated steel base, powder coated steel column.



#### Table Top

Table top is made of phenolic resin laminate with PVC cover.



#### Height Adjustment

Height adjustment by gas spring, used single-handed according to the convenience of the user.



#### Casters

Easy to install, more stable and firmer  
High quality medical caster, 1.5 inch /2 inch mute wheel.




#### Matching Design

It can be used with the bedside table, which is easy to disassemble. The combination of dining table and cabinet functions can serve two purposes and save space.





## CONSTANT INNOVATION




Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue